EU Lyme + VlsE IgG EIA Test System

EU Lyme + VlsE IgG Enzyme Immunoassay (EIA) Test System for the detection of IgG antibodies to *Borrelia afzelii*(PKO), *Borrelia garinii*(G2) and *Borrelia burgdorferi*(B31) including the VlsE Protein

Product # 44-8696G

96 Determinations

Store Kit at +2 to +8° C

INTENDED USE

The Trinity Biotech EU Lyme + VlsE IgG Enzyme Immunoassay (EIA) Test System is a qualitative test intended for use in the presumptive detection of human IgG antibodies to *Borrelia afzelii*, *Borrelia garinii*, and *Borrelia burgdorferi*(B-31) including the VlsE Protein in human serum. This EIA system should be used to test serum from patients with a history and symptoms of infection with *Borrelia*. All positive and equivocal specimens should be re-tested with a highly specific, second-tier test such as Western blot. Positive second-tier results are supportive evidence of infection with *Borrelia*. The diagnosis of Lyme disease should be made based on history and symptoms (such as erythema migrans), and other laboratory data, in addition to the presence of antibodies to *Borrelia afzelii*, *Borrelia garinii* and *B. burgdorferi*(B-31) including the VlsE protein. Negative results (either first or second-tier) should not be used to exclude Lyme disease.

I. SUMMARY

Lyme disease is a multisystem disease caused by the spirochete *Borrelia burgdorferi* (1). The disease has been documented in Europe since early this century. It was documented in the United States during an epidemic in 1975 among children in Old Lyme Connecticut who demonstrated arthritic symptoms. Steere, et al. recognized the disease as a unique clinical entity (2,3). The symptoms of Lyme disease have been mistaken with many diseases including: juvenile rheumatoid arthritis, lupus erythematosus, multiple sclerosis, Bell's palsy, rheumatic fever, Reiter's Syndrome, myocarditis and viral meningitis (4). The predominant strains present in Europe are *Borrelia afzelii*, *Borrelia garinii* and *Borrelia burgdorferi*.

The spirochete is transmitted by ticks of the genus *Ixodes* from animal reservoirs such as deer, mice, dogs, horses and birds. The ticks are commonly found on vegetation in endemic regions especially in wooded areas common to the animal reservoir. The incidence of human infection coincides with the tick season from May through September (3,5).

Although the symptoms of Lyme disease are varied and sometimes unclear, three distinct phases of the disease are recognized. Early manifestations include a single or multiple rash called erythema migrans (EM), a meningitis stage during the next weeks to months is often seen. Late manifestations are recognized to include arthritis or neurological signs and symptoms. In asymptomatic or subclinical cases, symptoms of infection may not be evident until the later stages of disease (5).

Isolation of *B. species* in culture is definitive evidence of active infection, but is not practical. Detection of specific antibodies is practical but an indirect marker of exposure. Patients produce IgM antibodies within a few weeks of the appearance of EM. Although only IgM antibodies may be detectable during the first month, IgG antibodies increase in most patients after approximately one month. Detectable levels of both IgG and IgM may persist for years (5,6).

Borrelia strains exhibit considerable antigenic variation. Patients often develop early antibodies to the flagellar antigen which can be cross reactive. Patients in the early stage of disease and a portion of patients with late manifestations may not have detectable antibodies. Early antimicrobial treatment, after appearance of EM may lead to diminished antibody concentrations. Serologic tests have been shown to have low sensitivity and specificity and, therefore, cannot be relied upon for establishing a diagnosis of Lyme disease (6,7,8).

The Second National Conference on Serological Diagnosis of Lyme Disease (1994) recommended the use of a two-tier test system for Lyme serology in which positive and equivocal samples from a sensitive first-tier test must be further tested by a more specific method such as Western blot (second tier). Positive results in the second tier test provide supportive evidence of exposure to *B. burgdorferi* which could support a clinical diagnosis of Lyme disease but should not be used as a criterion for diagnosis (9).

II. PRINCIPLE

The Trinity Biotech EU Lyme Disease IgG Test System is an indirect enzyme immunoassay (EIA) technique utilizing the purified native antigens of *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and *B. burgdorferi* (B 31) including the specific VlsE protein which are bound to polystyrene microwells. The antibody from the patient sample, which is added to the microwells in the first step of this procedure, binds with the antigen. Following a rinsing period which removes unbound antibody, a peroxidase labeled anti-human IgG is added to each microwell. The conjugate will attach during the second incubation step only if bound human IgG antibody is present from the first step reaction. Following a second rinsing period which removes unbound peroxidase conjugate, a color indicator solution is added to the microwells which will react only in the presence of bound peroxidase. An acid solution is added after a specified period of time in order to stop the enzymatic conversion of the indicator solution for spectrophotometric analysis.

III. REAGENTS SUPPLIED

Each kit contains the following components in sufficient quantities to perform the number of tests indicated on the package label.

1. Purified native antigens of *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and *B. burgdorferi* (B 31) antigens including the specific VlsE protein coated microassay plate: 96 wells, configured in a full plate, stored in a foil pouch with a desiccant. Product number 44-8208G
2. EU Lyme EIA IgG Positive Control: human serum. Proclin (0.1%) added as a preservative. Product number 44-8002, one vial, 100 µL.
3. EU Lyme EIA IgG Low Positive: human serum. Proclin (0.1%) added as a preservative. Product number 44-8003, one vial, 100 µL.
4. EU Lyme EIA Negative Control: human serum. Proclin (0.1%) added as a preservative. Product number 44-8001, one vial, 100 µL.
5. Peroxidase Conjugate, Anti-Human IgG. Goat anti-human IgG, containing Proclin (0.03%) as a preservative. Product number 44-8025G, one bottle, 13.5 ml.
6. EU Lyme EIA Color Developer: Tetramethylbenzidine (TMB). Product number 44-1006, one bottle, 13.0 ml.
7. Stop Solution: Contains a H₃PO₄ solution. Product number 40-1004, one bottle, 15 ml.
8. Serum Diluent. Contains Proclin (0.1%) as a preservative. Product Number 40-8012, one bottle, 100 ml.
9. 10x Wash Solution. Contains PBS, Tween-20. Product number 40-1013, one bottle, 100 ml.

IV. PRECAUTIONS

1. For *In Vitro Diagnostic* use.
2. The preservatives used in the reagents may be toxic if ingested. Do not mouth pipette.
3. This kit contains 1 molar phosphoric acid as a stop solution. Good laboratory procedure should be exercised when using acid solutions. Avoid contact with skin, eyes, or clothing. Use copious amount of water to wash exposed areas.
4. The human serum components used in the preparation of the Controls in this kit have been tested by an FDA approved method for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen and were found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
5. The Centers for Disease Control and Prevention, and the National Institutes of Health recommend that potentially infectious agents be handled at the Biosafety Level 2 (10).
6. Do not deviate from the specified temperature and timing requirements as listed in the package insert for both incubation and washing steps. Deviations may significantly alter the results of this test.
7. **All reagents must be brought to +20 to +25°C before performing this test procedure. Temperatures above or below the recommended range may result in substantial variation of the test results.**
8. If less than a full plate is to be used, center the strips within the plate frame to secure the strips during the wash steps.
9. Do not interchange kit components from one kit lot with another kit lot.
10. All unused microwells must be stored in desiccation at +2 to +8°C between uses.
11. Avoid cross-contamination of reagents. Wash hands before and after handling reagents. Cross-contamination of reagents and/or samples could cause false results.
12. Wells are to be washed three times.
13. If a sodium hypochlorite (bleach) solution is used as a disinfectant, do not expose the work area during the test procedure to avoid interference with the enzyme activity.
14. Assays designed for detection of antibodies to *Borrelia Afzelii*, *Borrelia Garinii* and *Borrelia burgdorferi* including the VlsE protein:
 - a. Do not indicate when exposure occurred.
 - b. Do not indicate that active replication is occurring or not occurring.
 - c. Frequently yield false-negative and false-positive results.
 - d. Do not have an established role in detecting such antibodies in serum collected during late manifestations of infection.
 - e. Do not often detect antibodies at clinical presentation during early Lyme disease, even when erythema migrans is present and isolation of *Borrelia species* in culture has been accomplished.
15. The concentrations of anti-*Borrelia* antibody in a given specimen determined by assays from different manufacturers may vary due to differences in assay methods and reagent specificity.
16. Do not mix components from different lots and do not interchange microwells or plates from IgG and IgM kits.

V. STORAGE AND STABILITY

1. Store kit at +2 to +8°C.
2. Bring all components to room temperature prior to use.
3. Refer to the expiration date on all reagents.
4. Refer to the expiration date on the microwell strips.
5. Immediately return kit to refrigeration after use.
6. All unused opened microwells must be stored in desiccation at +2 to +8°C between uses.

VI. ADDITIONAL MATERIALS REQUIRED BUT NOT SUPPLIED

1. Microtubes.
2. Microtube rack.
3. Storage container, 1.0 L, plastic or glass.
4. Pipetors, 10 µL, 100 µL, and 1.0 ml capacity.
5. Pipetor, 100 µL capacity, 8 channel (optional).
6. Waste container.
7. Timer, 30 minute.
8. Plastic squeeze bottle, 500 ml.
9. Absorbent paper (towels).
10. Distilled water.
11. EIA plate reader with 450 nm filter.
12. Desiccant storage for open microwell strips.

VII. REAGENT PREPARATION

1. Prepare 1X wash solution to be used in the assay from the 10X Wash Solution provided. On removal of the 10X Wash Solution from refrigeration, undissolved salts may be present. Allow the reagent to reach room temperature and shake the bottle to dissolve the salts. Do not dispense from the reagent until all salts are dissolved.
2. In order to prepare the 1X Wash solution, dilute 1 part of concentrate with 9 parts of distilled or deionized water in a clean plastic squeeze bottle. For each 8-well strip to be used, prepare a minimum of 50 ml diluted wash solution.
3. The Color Developer is ready to use.
4. The Conjugate is ready to use.
5. The Serum Diluent is ready to use.

VIII. SPECIMEN COLLECTION

Serological specimens should be collected under aseptic conditions. Avoid hemolysis by prompt separation of the serum from the clot. Serum should be stored at +2 to +8°C if it is to be analyzed within 4-7 days. Serum may be held for 3 to 6 months by storage at -20°C or lower. Lipemic and strongly hemolytic serum should be avoided. When specimens are shipped at ambient temperatures, additions of a preservative such as 0.01% thimerosal (merthiolate) or 0.1% sodium azide is strongly recommended. The NCCLS provides recommendations for storing blood specimens (Approved Standard Procedure for the Handling and Processing of Blood Specimens, H18-A2 1999) (11).

IX. PREPARATION OF SPECIMENS AND CONTROLS

A Positive Control, three Low Positives, a Negative Control, and a reagent blank must be run in each microplate tested. **Dilute kit controls at the same time as patient specimens for each run.** The reagent blank is prepared from Serum Diluent without the addition of human serum. Prepare the controls and patient specimens in appropriately marked pre-dilution tubes as follows:

1. Add 1.0 ml of Serum Diluent to each tube.
2. Add 10 µl of each serum to the appropriate tube and mix thoroughly.
3. Place each tube in an appropriate holding rack.
4. Specimens, controls, and calibrators are now diluted 1:100 and are ready for application to the *Borrelia* sensitized microwell plates.

X. TEST PROCEDURE

1. Remove the full plate from its foil pouch. If less than a full plate is used, strips must be centered within the plate frame to secure the strips during the wash steps. Place unused strips back into the foil pouch and seal tightly. Unused strips must be placed in desiccation at +2 to +8°C.
2. Add 100 µL of each control, specimen and reagent blank dilution into wells and record on EIA worksheet.
3. Incubate at (20-25°C) for 30 minutes.
4. Shake out contents of the wells. Refill each well with wash solution.
5. Repeat Step 4 two times. Blot strips dry on absorbent paper.
6. Add 100 µL of conjugate into each well.
7. Incubate at (20-25°C) for 30 minutes.
8. Repeat Steps 4 and 5.
9. Add 100 µL of Color Developer into each well.
10. Incubate at (20-25°C) for 15 minutes. Positive wells will develop a blue color.
11. Add 100 µL of Stop Solution to each well. Mix by gently tapping side of plate to evenly disperse Stop Solution. Do not create bubbles. Positive wells will develop a yellow color.
12. Incubate for 2 minutes.
13. Read the test strips by an EIA plate reader at 450nm within 30 minutes after addition of stop solution.

XI. QUALITY CONTROL

A spectrophotometer or EIA plate reader with a 450 nm filter is needed for the analysis of the color intensity.

1. Read the plate at 450 nm. Print out the optical density values (OD₄₅₀).
2. Round the blank well OD to 2 digits to the right of the decimal and subtract from all test ODs. (0.795 = 0.80, 2.333 = 2.33).
3. Round all test ODs as above.
4. Calculate the mean OD of the three Low Positive wells. If any of the three wells has an OD blank more than 10% from the mean, disregard that value and recalculate the mean from the remaining two wells. If the ODs from the remaining 2 wells are not within 10% of the mean, assay results are invalid and must be repeated.
5. To calculate a Lyme Index Value (LIV), divide the OD of the sample by the mean O.D. of the Low Positive Control.
6. The assay is considered valid and reportable when:
 - a. The Positive Control has a Lyme Index Value ≥ 1.50 .
 - b. The Negative Control has a Lyme Index Value < 0.50 .
 - c. The Low Positive Control has an O.D. value of 0.45-0.95.
 - d. The Blank has an O.D. value of < 0.175 .

XII. INTERPRETATION OF RESULTS

A Lyme Index Value of 1.0 or greater should be considered presumptive evidence of exposure to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and/or *B. burgdorferi* (B 31) including the specific VlsE protein

XIII. REPORTING**Positive - Lyme Index of ≥ 1.0**

A positive result is presumptive evidence of exposure to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and/or *B. burgdorferi* (B 31) including the specific VlsE protein. The result should be supplemented by a second tier test that is more specific for antibodies to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and *B. burgdorferi* (B 31) including the specific VlsE protein and should not be reported until second-step testing is complete. A positive result is not an indication of when exposure occurred.

Negative - Lyme Index of < 0.5

No detectable antibodies to *B. burgdorferi*. A negative result indicates that there was not serologic evidence of exposure to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and/or *B. burgdorferi* (B 31) including the specific VlsE protein at the time the specimen was collected. A negative result should not be the basis for excluding *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and *B. burgdorferi* (B 31) including the specific VlsE protein as the cause of illness, especially if blood was collected within 2 weeks of when symptoms began. If Lyme disease is strongly suspected, a second specimen should be collected 2 to 4 weeks after the first specimen and tested then.

Equivocal - Lyme Index of ≥ 0.5 and ≤ 0.99

An equivocal result is presumptive evidence of possible exposure to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and/or *B. burgdorferi* (B 31) including the specific VlsE protein. The result should be supplemented by a second tier test that is more specific for antibodies to *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and *B. burgdorferi* (B 31) including the specific VlsE protein and should not be reported until second-step testing is complete.

XIV. LIMITATIONS OF PROCEDURE

1. Sera from patients with other pathogenic spirochetal diseases such as syphilis, yaws, pinta, leptospirosis, and relapsing fever may give false positive results. Sera from patients with mononucleosis or lupus erythematosus may also give false positive results. In cases where false positive results occur, clinical, epidemiologic and laboratory workups should be carried out. Although the clinical picture is quite different between active syphilis and Lyme disease, an easy means of differentiating these two diseases is by the use of the VDRL or RPR tests. In active syphilis, the VDRL or RPR are positive. In Lyme disease, the VDRL and RPR are negative.
2. Patients with early stages of Lyme disease may not test positive with this test. The IgG antibodies may not have reached detectable levels in early Lyme during the EM clinical phase of the disease. Negative results in early stages of disease have low predictive value.
3. In order to test for early stage I infection, use a reagent system designed to detect IgM antibody without interference from rheumatoid factor and IgG antibodies.
4. Antibiotic therapy given early in the disease may prevent the development of an antibody response.
5. The evaluation of all test results must include the clinical history presented by the patient, the patient's exposure in endemic regions for Lyme disease, epidemiological data, other test results, and any other spirochetal diseases. Positive test results in patients for which a history of exposure to Lyme disease, or symptoms or clinical findings consistent with Lyme disease do not exist may have a low predictive value.
6. Positive or equivocal first-tier test results should not be reported until second-tier testing of the specimen is performed using a method that is more specific such as Western blot.
7. The use of this assay has not been evaluated for individuals who have received Lyme disease vaccine.

XV. EXPECTED VALUES

In patients without present or past infection with *Borrelia afzelii* (PKO), *Borrelia garinii* (G2) and/or *B. burgdorferi* (B 31), test results should be negative except for patients with cross-reacting antibodies (See Limitations).

In patients with Lyme disease, the test results are dependent on the stage of the disease. IgG antibody titers can be detected within two or three weeks following onset of disease. The VlsE protein is a very early marker of the IgG immune response and antibodies to the VlsE may persist. As the disease progresses, the chances of a positive test increases. Persons presenting within 3 weeks of onset of EM are positive for IgG approximately 25% of the time. Persons with second stage symptoms are almost always positive as are persons with third stage symptoms.

IgG testing is most useful for the detection of ongoing or past cases of Lyme disease. For early detection, a test which detects the IgG response to the VlsE protein or IgM class antibodies is recommended. IgM antibodies rise quickly after infection, reaching their peak within 3 weeks and declining rapidly. IgG antibodies start rising 2 to 3 weeks after infection and persist while symptoms are present and drop slowly during recovery (5,6).

XVI. SPECIFIC PERFORMANCE CHARACTERISTICS**Specificity**

The Specificity of Trinity Biotech IgG was evaluated at an independent diagnostic laboratory in England. The study was conducted on a panel of 100 serum samples taken from presumably healthy donors living in an area endemic for Lyme borreliosis.

Relative Specificity 98% (98/100)

Sensitivity

The study was conducted on three panels of serum samples from patients with some of the most common clinical manifestations of a *Borrelia* infection: 55 serum samples from patients with Erythema migrans, 30 serum samples from patients with Arthritis and 40 serum samples from patients with Neuro-Borreliosis.

	Diagnostic Sensitivity for IgG	Diagnostic Sensitivity between both IgG and IgM
Erythema migrans	9 out of 55 84%	8 out of 55 85.5%
Arthritis	3 out of 30 90%	2 out of 30 93%
Neuro-Borreliosis	1 out of 40 97.5%	0 out of 40 100%

Cross-Reactivity

Sera from patients with syphilis, inflammatory diseases (rheumatoid factor (RF) positive, autoimmune and *Toxoplasma gondii*) were tested in Trinity Biotech IgG with the following results:

Patient Sera	Number of sera	Number of Positives
Syphilis	20	2
RF	20	1
Autoimmune	20	0
Toxoplasma gondii	20	0

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MarDx Diagnostics, Inc.