ISLET CELL ANTIBODY (ICA_b) TEST SYSTEM

IVD

PRODUCT INSERT

REF	1123 ICA (Primate Pancreas) 40 Determinations
-----	---

INTENDED USE

An indirect immunofluorescence test for the detection and semi-quantitation of islet cell antibodies in human serum to aid in the diagnosis of *insulin dependent diabetes mellitus* (IDDM).

SUMMARY AND EXPLANATION

Diabetes is a chronic and complex metabolic disease influenced by various hereditary and environmental factors that result in the inability of the body to maintain the use of carbohydrates, fats and proteins. The condition, characterized by high blood glucose levels, is caused by a deficiency in insulin production or an impairment of insulin utilization. Most cases of diabetes fall into two clinical categories: *insulin-dependent diabetes mellitus* (IDDM or Type I diabetes) and *non-insulin dependent diabetes mellitus* (NIDDM or Type II diabetes). Prognosis, treatment and disease management are different for each type. It is well accepted that IDDM is an autoimmune disease targeting β -cells of the islets of Langerhans in the pancreas. The autoimmune response to islet cell antigens elicits antibody responses to antigens such as glutamic acid decarboxylase (GAD), ICA-512 and insulin. They have been found to be highly predictive markers, particularly if present in high titer¹⁻¹⁵. Detection of these ICA_b's by indirect immunofluorescence (IF) on pancreas substrate is considered the gold standard for diagnosis of IDDM^{16,17}. These cytoplasmic ICA_b's are currently used for the prediction of type I diabetes¹⁸⁻²⁰. ICA_b's are detected in up to 90% of newly diagnosed diabetic patients. In the Bart's-Windsor family study, 100% of the first degree relatives of IDDM patients with ICA_b's >80 JDF units progressed to IDDM within 10 years. The level of ICA_b's appears to be highest prior to the onset of Type I diabetes and diminishes progressively thereafter^{12,18}. ICA_b's have several distinct specificities and show two major patterns of reactivity⁴. The first pattern is restricted predominantly to the β -cells. The second stains all cells within the islet, and is the classical staining pattern for cytoplasmic ICA_b's.

PRINCIPLES OF PROCEDURES

Islet cell antibodies are detected by indirect immunofluorescence (IF). In this method, patients' sera are incubated on sections of monkey pancreas to allow binding of antibodies to the tissue substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of these sections with fluorescein-labeled conjugates. Reactions are observed under a fluorescence microscope equipped with appropriate filters.

The presence of islet cell antibodies is demonstrated by an apple green fluorescence of the cytoplasm of the islets of Langerhans. The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions of the positive control²¹.

PRODUCT INFORMATION

Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

Materials Provided

ImmuGlo™ ICA_b Primate Pancreas

REF	1123
-----	------

Kits contain sufficient reagents to perform 40 determinations each.

10 x

SORB	SLD	4
------	-----	---

4 well **Monkey Pancreas Substrate Slides**

1 x 0,5 ml

CONTROL	+	ICA
---------	---	-----

ICA Positive Control. Human serum containing islet cell antibodies, standardized against the JDF (Juvenile Diabetes Foundation) positive reference serum. **JDF units are stated on the vial label.**

EN

1 x 0,5 ml **CONTROL -**

1 x 5 ml **IgG-CONJ FITC EB**

1 x 5 ml **CONJ B**

1 x 60 ml **BUF ICA**

2 vials **BUF WASH**

1 x 5 ml **MOUNTING MEDIUM**

1 x 12 ml **COVER SLD**

Optional Components

1 x 5 ml **IgG-CONJ FITC**

1 x 1 ml **EVANS**

1 x 5 ml **IgG-CONJ FITC PA**

Negative Control. Contains human serum.

Anti-human IgG FITC Conjugate containing Evan's Blue. Protect from light.

Conjugate B. Protect from light.

ICA Sample Diluent.

Phosphate Buffered Saline (PBS). Dissolve each vial to 1 liter.

Mounting Medium. Do not freeze.

Coverslips.

REF 2100X - Anti-human FITC Conjugate. Protect from light.

REF 2510* - Evan's Blue Counterstain.

REF 2099 - Anti-human IgG FITC Conjugate. Protect from light. Conjugate minimizes background reactions on primate tissues.

Symbols used on labels:

LOT Lot number

REF Catalog number

 Use by

 Storage temperature

 Consult instructions for use

IVD In vitro diagnostic use

 Manufacturer

 Number of Tests

 Date of Manufacture

Danger. May cause cancer. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. IF exposed or concerned: Get medical advice/attention. Wear protective gloves, protective clothing, eye protection, face protection. Store locked up. Dispose of contents/container to comply with local, state and federal regulations.

Materials Required But Not Provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation Chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials²².

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident

immediately to laboratory director or poison control center. Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Do not use beyond expiration date.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum 1:5 with the Sample Diluent provided (0.1 ml serum + 0.4 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (approximately 50 µl of the Negative Control to well #1. Similarly apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides 18 - 24 hours at 2°-8°C.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the IgG FITC Conjugate dropper vial and gently squeeze to apply 1 drop (approximately 50 µl) to each well.
9. Place the lid on the incubation chamber and incubate slides 30 minutes at room temperature.
10. Repeat steps **7 through 9** for each slide, using Conjugate B in place of IgG FITC Conjugate.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with PBS for 10 minutes. If optional conjugate without counterstain is used (see optional components in Materials Provided Section), 2-3 drops of Evan's Blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.

EN

12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Each test run should include the Positive and Negative Controls. Make serial two-fold dilutions starting at 1:5. The reciprocal of the highest dilution producing a positive reaction is the titer. If the Positive Control titer is within the limits defined by the enclosed QC specifications, the levels of the antibody in patient serum can then be reported in JDF units¹⁶. The JDF unit value of the Positive Control is printed on the vial label. To calculate the JDF unit value of the unknown serum, simply divide the titer of the unknown by the titer of the Positive Control and multiply this by the JDF units on the Positive Control label.

Example:

Positive Control titer is 1:4. Unknown Sample titer 1:10. Positive Control label reads 160 JDF units.

Calculation:

$$\text{ICAb Concentration} : \frac{10}{4} \times 160 = 400 \text{ JDF Units}$$

Preparation of Serial Dilutions

Number four tubes 1 through 4. Add 0.4 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 4. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4
Serum	0.1 ml			
	+			
Buffered Diluent	0.4 ml	0,2 ml	0,2 ml	0,2 ml
		↗	↗	↗
Transfer		0,2 ml	0,2 ml	0,2 ml
Final dilution	1:5	1:10	1:20	1:40 etc.

QUALITY CONTROL

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence of the of the islet cells, whereas the Positive Control should stain the cytoplasm of the islet cells. If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include : improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

EN

RESULTS

The results of the tests for islet cell antibodies should be reported as negative (<5) or positive with titer or JDF units⁸. Read for specific staining of the cytoplasm of the isletofLangerhans. Various other tissue antibodies such as antinuclear antibodies (ANA), mitochondrial antibodies, and smooth muscle antibodies may also be observed on pancreas sections. Sera giving any of these reactions should be reported as negative for antibodies to islet cells. Any sera giving nuclear staining reactions may be tested on HEp-2 cells or mouse liver sections. Any sera giving smooth muscle or mitochondrial staining reactions may be tested on mouse kidney/ stomach sections. See ICAb positive reaction in Photo 1 at the end of this document.

NOTE: ICAb reactions are, by their nature, much weaker than reactions for ANA or most other immunofluorescence antibody reactions.

LIMITATIONS OF THE PROCEDURE

Occasionally sera may exhibit strong staining for ANA or other autoantibodies. These may interfere with the ability to detect islet cell antibodies. In such cases, titrating the serum may permit the visualization of the islet cell antibodies. In other cases their titer may be lower than the ANA or other antibodies and hence may not be detected. The ImmuGlo™ islet cell antibody test should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only.

EXPECTED VALUES

Approximately 50-80% of new-onset diabetics are positive for islet cell antibodies. The prevalence of islet cell antibodies in non-diabetic first-degree relatives and in non-diabetic⁵⁻¹⁶ normal subjects has been reported to be 2-5% and 0.25-1.7%, respectively. Approximately 11% of islet cell antibody positive first-degree relatives develop diabetes each year. Intervals of as long as 8 years between detection of islet cell antibodies and the onset of diabetes have been reported. See ICAb incidence in Table 1 at the end of this document.

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REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA

1. Riley WJ, Maclaren NK, Krischer J et al. N. A prospective study of the development of diabetes in relatives of patients with insulin dependent diabetes. *New Engl J Med*; 1990, 323:1167-1172.
2. Neufeld M, Maclaren NK, Riley WJ et al. Islet cell and other organ-specific antibodies in U.S. Caucasians and Blacks with insulin-dependent diabetes mellitus. *Diabetes*; 1980, 29:589-592.
3. Gianani R, Pugliese A, Bonner-Weir S et al. Prognostically significant heterogeneity of cytoplasmic islet cell antibodies in relatives of patients with type I diabetes. *Diabetes*; 1992, 41:347-353.
4. Ziegler AG, Herskowitz RD, Jackson RA et al. Predicting Type I Diabetes. *Diabetes Care*; 1990, 13:762-775.
5. Gale EAM and Bottazzo GF. Can we predict type I (insulin dependent) diabetes? In "World Book of Diabetes in Practice"; 1986, Vol 2, Krall L, Ed, Elsevier, New York, 25-29.
6. Gorsuch AN, Spencer KM, Lister J et al. Evidence for a long pre-diabetic period in Type I (insulin dependent) diabetes mellitus. *Lancet*; 1981, 2:1363-1365.
7. Tarn AC, Bonifacio E, Dean BM et al. Predicting insulin-dependent diabetes (Letter). *Lancet*; 1988, 2:627-628.
8. Jackson RA, Soeldner JS and Eisenbarth GS. Predicting insulin-dependent diabetes (Letter). *Lancet*; 1988, 2:627-628.
9. Srikanta S, Ganda OP, Gleason RE et al. Pre-type I diabetes: linear loss of beta cell response to intravenous glucose. *Diabetes*; 1984, 33:717-720.
10. Srikanta S, Ganda OP, Rabizadeh A et al. First degree relatives of patients with type I diabetes mellitus: islet cell antibodies and abnormal insulin secretion. *N Engl J Med*; 1985, 313:461-464.
11. Srikanta S, Ganda OP, Jackson RA et al. Pre-type I diabetes: common endocrinologic course despite immunologic and immunogenetic heterogeneity. *Diabetologia*; 1984, 27:146-149.
12. Riley WJ, Spillar RP, Waltz J and Brody B. Predictive value of islet cell antibodies (ICA) - 6 years experience (Abstract). *Diabetes*, 1983, 33 (Supp 1):44A.
13. Riley W, Maclaren N, Spillar RP et al. Predictive value of ICA for IDD and insulinopenia to iv glucose (Abstract): *Diabetes* 37 1988, (Issue 5 Suppl): 5A.
14. Maclaren NK, Horne G, Spillar RP et al. Islet cell antibodies (ICA) in U.S. school children (Abstract). *Diabetes*; 1985, 34 (Suppl 1):84A.
15. Spencer KM, Tarn A, Dean BM et al. Fluctuating islet cell autoimmunity in unaffected relatives of patients with insulin dependent diabetes. *Lancet*; 1984, 1:764- 766.
16. Greenbaum J, Palmer JP, Nagataki S et al. Improved specificity of ICA assays in the fourth international immunology of diabetes serum exchange workshop. *Diabetes*; 1992, 41:1570-1574.
17. Bonifacio E, Lemmark Å, Dawkins RL et al. Serum exchange and use of dilutions have improved precision of measurement of islet cell antibodies. *J Immunol Methods*; 1988, 106:83-88.
18. Kolb H, Dannehl K, Grueneklee D et al. Prospective analysis of islet cell antibodies in children with type I (insulin dependent) diabetes; 1988, *Diabetologia*; 1988, 31:189-194.
19. Srikanta S, Ganda OP, Jackson RA et al. Type I diabetes mellitus in monozygotic twins: chronic progressive cell destruction. *Ann Intern Med*; 1983, 99:320-326.
20. Betterle C, Presotto F, Pedini B et al. Islet cell and insulin autoantibodies in organ specific autoimmune patients: their behavior and predictive value for the development of type I diabetes mellitus: a 10 year follow-up study. *Diabetologia*; 1987, 30:292-297.
21. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York; 1987, 3rd Ed, 3-40.
22. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health. [HHS Pub. No. (CDC) 1999, 93-8395].

Photo 1. ICAb Staining Reaction

Note specific staining of the cytoplasm of the islet of Langerhans.


Table 1. Incidence of Anti-Islet Cell Antibodies


Disease Group	Age (years)	No. Patients examined	% Positive	
Type I Diabetes (IDDM)	at onset	<1-10	19	63
		11-20	25	60
		21-40	8	25
	long standing	<1-10	22	41
		11-20	71	39
		21-40	26	24
		41-70	13	0
	71-80	3	33	
Type II Diabetes (NIDDM)	at onset	<1-40	0	-
		41-80	39	3
	long standing	<1-10	0	-
		11-20	5	20
		21-80	75	1
	Non-diabetic first-degree relatives	<1-30	61	0
31-50		119	2	
51-80		19	0	
non-diabetic controls	>18	200	0	

For technical assistance please contact:

Immco Diagnostics, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands