## Anti-Glomerular Basement Membrane (GBM) Test System

**IVD**

### PRODUCT INSERT

**REF** 1124 Anti-GBM (Primate Kidney) 48 Determinations

#### INTENDED USE

An indirect immunofluorescence (IFA) antibody test for the detection and semi-quantitation of anti-glomerular basement membrane (GBM) antibodies in human serum.

#### SUMMARY AND EXPLANATION

Rapidly progressive glomerulonephritis (RPGN) is a clinical syndrome developing over days or weeks characterized by crescentic glomerulonephritis as seen by histopathology of the kidney. The prognosis is poor if not recognized early and if an appropriate treatment is not instituted. RPGN may be assessed based on indirect immunofluorescence serum studies and electron microscope evaluations of renal biopsies.

Using the above criteria RPGN may be classified into:

- a. Immune complex mediated disease characterized by the presence of anti-DNA antibodies or anti-streptococcal antibodies.
- b. Anti-glomerular basement membrane (GBM) mediated glomerulonephritis and Goodpasture's syndrome.
- c. Anti-neutrophil cytoplasmic antibody (ANCA) associated glomerulonephritis.

In a study by Jayne et al<sup>1</sup>, of 889 RPGN suspected patients, 47 (5%) had anti-GBM, 246 (28%) had ANCA and 576 (65%) had neither antibodies. 2% had both ANCA and anti-GBM antibodies. The frequency of anti-GBM antibodies and associated diseases are listed in Table 1.

#### PRINCIPLES OF PROCEDURES

In this IFA method, to enhance the specific fluorescence staining, the tissue substrate is first incubated with an Antigen Enhancing Buffer prior to the incubation with patient sera. Patients' sera are incubated on monkey kidney substrate to allow binding of antibodies to the substrate. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgG class are detected by incubation of the substrate with fluorescein-labeled, anti-human IgG conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters.

The presence of anti-GBM antibody reactions are demonstrated by an apple green fluorescence of the basement membrane of the kidney glomeruli. The titer (the reciprocal of the highest dilution giving a positive reaction) is then determined by testing serial dilutions.

#### PRODUCT INFORMATION

##### Storage and Preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

##### Materials Provided

ImmuGlo™ Anti-GBM (Primate Kidney) **REF** 1124

Kits contain sufficient reagents to perform 48 determinations each.

## EN


8 x	<b>SORB</b> <b>SLD</b> <b>6</b>	6 well <b>Monkey Kidney Substrate Slides</b>
1 x 0.5 ml	<b>CONTROL</b> <b>+</b> <b>GBM</b>	<b>GBM Positive Control.</b> Human serum containing GBM antibodies.
1 x 0.5 ml	<b>CONTROL</b> <b>-</b> <b>GBM</b>	<b>GBM Negative Control.</b> Contains human serum.
1 x 5 ml	<b>IgG-CONJ</b> <b>FITC</b> <b>EB</b> †	Anti-human IgG FITC Conjugate containing Evan's Blue. <b>Protect from light.</b>
1 x 5 ml	<b>DIL-ENH</b> <b>GBM</b>	<b>Antigen enhancing buffer.</b>
1 x 60 ml	<b>BUF</b> <b>GBM</b>	<b>GBM sample diluent.</b>
2 x	<b>BUF</b> <b>WASH</b>	<b>Phosphate Buffered Saline (PBS).</b> Dissolve each vial to 1 liter.
1 x 5 ml	<b>MOUNTING</b> <b>MEDIUM</b>	<b>Mounting Medium.</b> Do not freeze.
1 x 12	<b>COVER</b> <b>SLD</b>	<b>Coverslips.</b>

### Optional Components

1 x 5 ml	<b>IgG-CONJ</b> <b>FITC</b> †	<b>Anti-human FITC Conjugate. Protect from light.</b>
1 x 1 ml	<b>EVANS</b>	<b>Evan's Blue Counterstain*</b>

† Replaces conjugate without counterstain in Code numbers containing EB

Symbols used on labels:

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Consult instructions for use
	Number of tests
	Manufacturer
	Date of Manufacture

 \* Danger. May cause cancer. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. IF exposed or concerned: Get medical advice/attention. Wear protective gloves, protective clothing, eye protection, face protection. Store locked up. Dispose of contents/container to comply with local, state and federal regulations.

### Materials Required But Not Provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes
- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

**WARNINGS AND PRECAUTIONS**

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials<sup>2</sup>.

**WARNING** – Sodium azide (NaN<sub>3</sub>) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

Instructions should be followed exactly as they appear in this insert to ensure valid results. Do not interchange kit components with those from other sources other than the same catalog number from Immco Diagnostics Inc. Do not use beyond expiration date.

**SPECIMEN COLLECTION AND HANDLING**

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum should be frozen at -20°C. Avoid repeated freezing and thawing of samples.

**PROCEDURE****Test Method****A. Screening**

1. Dilute each patient serum **1:10** with the Buffered Diluent provided (0.1 ml serum + 0.9 ml Diluent). Do not dilute Positive or Negative Controls. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for **10-15 minutes**. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Pipet **1 to 2 drops** of **Antigen Enhancing Buffer** onto each well. Place the lid on the incubation chamber and incubate **30 minutes** at room temperature.
5. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash **10 minutes**. Repeat process with all remaining slides.
6. Remove slide from Coplin jar and shake off excess liquid and place the slide in the incubation chamber. Gently place blotter onto slide and press.
7. Invert dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) of the **Negative Control** to well #1. Similarly apply **1 drop** of **Positive Control** to well #2. Avoid overfilling the wells.
8. Using a micropipette or Pasteur pipette, apply **1 drop** of patient's diluted serum (approximately 50 µl) to the other wells. Avoid overfilling the wells.
9. Place the lid on the incubation chamber and incubate **30 minutes** at room temperature.
10. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash **10 minutes**. Repeat process with all remaining slides.
11. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the **Conjugate** dropper vial and gently squeeze to apply **1 drop** (approximately 50 µl) to each well. Repeat process with all remaining slides.
12. Replace the lid on the incubation chamber and incubate **30 minutes** at room temperature.
13. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately **10 ml** of PBS using a pipette, or rinse slide in a beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash **10 minutes**. Repeat process with all remaining slides. **NOTE:** Improper washing may lead to increased background fluorescence.
14. Remove a slide from the Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. **While slide is still wet mount the coverslip.** Place **3 drops** of the Mounting Medium evenly spaced on a coverslip and invert the slide onto the coverslip. To remove any air bubbles gently apply pressure along the edge of the coverslip. Avoid any movement of the coverslip. Repeat process with all remaining slides.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of **200x** or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8°C.

**B: End Point Determination (Titration)**

A serial two-fold dilutions starting at 1:10 (see below). Using one slide, a serum may be tested at dilutions ranging from 1:10 to 1:320. If positive at a 1:320 dilution, the titer is reported as greater or equal to 320. Additional slides may be used to obtain

## EN

endpoints for those sera still positive at a 1:320 dilution. The reciprocal of the highest dilution producing a positive reaction is the titer.

### Preparation of Serial Dilutions

Number six tubes 1 through 6. Add 0.9 ml of Buffered Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

<b>Tubes</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>
<b>Serum</b>	0.1 ml					
	+					
<b>Buffered Diluent</b>	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↻	↻	↻	↻	↻
<b>Transfer</b>	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	
<b>Final dilution</b>	1:10	1:20	1:40	1:80	1:160	1:320 etc.

### QUALITY CONTROL

Both Positive and Negative Controls should be included with each test run. The negative control should show no specific fluorescence of the kidney glomeruli, whereas the positive control should have 2+ or greater staining intensity of the basement membrane of the kidney glomeruli.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control.
- Problems with the optical system of the fluorescence microscope. These may include: improper alignment, use of the bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

### RESULTS

The results of the tests for anti-GBM should be reported as negative (<10), positive greater or equal to 320, or preferably, positive with titer. Read for specific linear basement membrane staining of the kidney glomeruli. See photo 1 at the end of this document.

### LIMITATIONS OF THE PROCEDURE

In some cases, sera positive for anti-GBM antibodies may either be very weak or negative at the initial screening dilution (prozone phenomenon). In such doubtful cases the sera should be screened at higher dilutions and, if positive, antibody titers determined.

In some cases the presence of two or more antibodies in a serum which are reactive with the same substrate may cause an interference in their detection by immunofluorescence. This interference may cause either a failure to detect anti-GBM antibodies or a suppression of the titer if the interfering antibody has a higher titer than anti-GBM antibodies.

### EXPECTED VALUES

Anti-GBM antibodies react with the glomerular basement membrane, the tubular basement membrane or the alveolar basement membrane with a continuous and linear pattern<sup>3</sup>. The anti-GBM antibodies are primarily of the IgG class. Approximately 5% of the patients with rapidly progressive glomerulonephritis (RPGN) are anti-GBM antibody positive<sup>1</sup>. Goodpasture's Syndrome (glomerulonephritis and pulmonary hemorrhage) accounts for two thirds of the patients positive for anti-GBM antibodies. RPGN and milder forms of glomerulonephritis account for the other one third. Occasionally, patients with anti-GBM antibodies may have only lung involvement<sup>4</sup>.

Anti-GBM antibodies are associated with Goodpasture's syndrome and RPGN. Circulating anti-GBM antibodies are present in 90% of patients with Goodpasture's syndrome and 60% of patients with RPGN. These antibodies may disappear 4-8 months after bilateral nephrectomy<sup>5</sup>.

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
**REFERENCES • ΒΙΒΛΙΟΓΡΑΦΙΑ • LITERATUR • BIBLIOGRAPHIE • BIBLIOGRAFIA**

1. Jayne DRW, Marshall PD, Jones SJ and Lockwood CM. Autoantibodies to GBM and neutrophil cytoplasm in rapidly progressive glomerulonephritis. *Kidney Int* 37:965-970, 1990.
2. Biosafety in Microbiological and Biomedical Laboratories. Center for Disease Control, National Institute for Health. (HHS Pub. No {CDC} 93-8395) 1999.
3. Wilson CB. Radioimmunoassay for anti-glomerular basement membrane antibodies. In "Manual of Clinical Immunology," Rose N and Friedman H (Eds), 2nd ed, ASM, pp 376-379, 1980.
4. Wilson CB and Dixon FJ. Renal injury from immune reactions involving antigens in or of the kidney. In "Contemporary Issues in Nephrology", Wilson CB, Brenner BM and Stein JH (Eds), Vol 3, Churchill Livingstone, New York, pp 35-66, 1979.
5. Jenis EM and Lowenthal DT. *Kidney Biopsy Interpretation*. FA Davis Co. Philadelphia, 1977.
6. Fish AJ, Kleppel M, Jeraj K and Michael AF. Enzyme immunoassay of anti-glomerular basement membrane antibodies. *J Clin Med* 105:700-705, 1985. 11

**Photo 1. Positive GBM staining reaction on monkey kidney, 200X. Note staining of the basement membrane of the glomeruli.**

**Table 1. Frequency of Anti-GBM Antibodies**

<b>Disease Group</b>	<b>% Positive</b>
Anti-GBM nephritis	100
Anti-TBM nephritis	0
Vasculitis	0
Anaphylactoid purpura nephritis	0
Systemic lupus erythematosus	0
Membranoproliferative glomerulonephritis	0
IgA nephropathy	0
Normals	0

Adapted from Fish AJ et al.<sup>6</sup>

*For technical assistance please contact:*

**IMMCO Diagnostics, Inc.**  
**60 Pineview Drive**  
**Buffalo, NY 14228-2120 USA**  
**Telephone: (716) 691-0091**  
**Fax: (716) 691-0466**  
**Toll Free USA/Canada: 1-800-537-TEST**  
**E-Mail: [info@immco.com](mailto:info@immco.com)**

*or your local product distributor*

**EC|REP**

**EMERGO EUROPE**  
**Prinsessegracht 20**  
**2514 AP The Hague**  
**The Netherlands**