ImmuBlot™

Anti-P0 Antibody

Western Blot

IVD

PRODUCT INSERT

REF 1192

20 Determinations

INTENDED USE

A Western Blot immunoassay for use in the detection and identification of antibodies to P0, a 30kD antigen associated with autoimmune hearing loss.

SUMMARY AND EXPLANATION

Hearing loss may be caused by a number of conditions. Certain types of hearing loss can be reversed if diagnosed early and appropriate treatment is instituted. Sensorineural hearing loss (SNHL), commonly referred to as nerve deafness, may be caused by genetic or acquired factors such as infections or can be immunologically mediated. Correct diagnosis of SNHL can be made with the aid of a combination of comprehensive patient history and laboratory studies. In the majority of cases, no cause of SNHL is apparent: such cases are referred to as idiopathic SNHL. Subgroups of idiopathic SNHL cases are treatable with immunosuppressive therapy with gratifying results^{1, 2}. The laboratory work up to identify these cases may include serum antibody tests to 68 KD (hsp-70) inner ear antigen³⁻⁸ and other auto antigens such as to P0, a 30 KD antigen⁷⁻⁹. Approximately 60% of patients with progressive SNHL have antibodies to P0. A strong positive correlation between chronic progressive hearing loss and the presence of anti-P0 antibodies was observed. In addition patients with positive P0 antibodies had bilateral hearing loss.

PRINCIPLE OF PROCEDURE

This Western Blot immunoassay uses purified P0 antigen from the bovine or primate tissue. To perform the test, PVDF strips are incubated with diluted patient serum. Antibodies specifically bind to the P0 protein on the strip. After proper washing and incubation steps with the Conjugate, strips are washed and incubated with a precipitable substrate. Anti-P0 antibody positive reactions appear as a blue-violet band at 30kD, sometime this reaction is seen as a doublet.

REAGENTS

Storage and Preparation

Store all reagents at 2-8°C. **Do not freeze.** Do not use if liquid reagents are turbid or a precipitate is present. Prior to starting the assay, reagents must be equilibrated to room temperature (~22°). Antigen strips can only be used once. Do not interchange components of different lots. Do not use reagents beyond expiration date indicated on labels.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by US FDA required tests. However, all human blood derivatives and patient specimens should be considered potentially infectious and good laboratory practices in storing; dispensing and disposing of these materials must be followed¹⁰.

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. NaN₃ is toxic if ingested. Report incidents immediately to laboratory director or poison control center.

Follow good laboratory practices to minimize microbial and cross contamination of reagents.

Materials provided**ImmuBlot™ Anti-P0 Antibody 1192**

The ImmuBlot™ anti-P0 (30kD) antibody Western Blot kit*, contains sufficient reagents to perform 20 determinations.

1 x 18	Western Blot Strips, with P0 antigen.
1 x 120 µl	Positive Control (<i>red vial cap</i>)
1 x 120 µl	Negative Control (<i>yellow vial cap</i>).
1	Control Card
1 x 250 µl	Conjugate (<i>blue vial cap</i>)
1 x 60 ml	Serum Diluent
1 x 25 ml	Enzyme Substrate (<i>amber bottle</i>)
1 vial	Powdered Wash Buffer; reconstitute to one liter with deionized or distilled water.
3	Assay trays

Symbols used on labels

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Consult instructions for use
	Number of tests
	Manufacturer
	Date of Manufacture

Materials Required But Not Provided

- Clean 1000 ml graduated cylinder
- Non-serrated forceps (Filter forceps)
- Rocker or rotating platform shaker
- Absorbent paper or paper towels
- Deionized or distilled water
- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 10 to 1000 µl
- Disposable pipette tips

- Timer

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbial contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2-8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes

- Read Product Insert carefully before starting with the assay.
- Let serum specimens and test reagents equilibrate to room temperature for ~30 minutes prior to starting the test procedure. Return all unused specimens and reagents to the refrigerator promptly after use.
- Proper washing technique is critical to the satisfactory performance of the assay.
- Manipulate test strips with clean forceps only. Do not touch with bare hands.
- Strips are individually numbered at the top of each strip. Assign specimen identification numbers to the respective strips on the Report Form.
- Complete all other relevant information on the Report Form prior to starting the assay.

Test Method

- Step 1** Using blunt forceps, place required number of **Strips** labeled side up into individual wells of the assay tray.
- Step 2** Pipette **1.0 ml** of Serum Diluent into each well.
- Step 3** Pipette **10 µl** of Positive and Negative Control and patient sample into appropriate wells to obtain a **1:101 dilution**. Incubate **60 minutes** (±5 min.) at room temperature on a rocker or rotating shaker.
- Step 4** Aspirate sample solution into waste container. Thoroughly wash strips with Wash Buffer by squirting approximately 2ml of solution directly onto strips. Wash strips with gentle agitation for **5 minutes** and aspirate solution into waste container. **Repeat 2x**. *Caution: Complete washing of the strips between incubations is crucial to obtain valid results. Improper washing will result in high background staining.*
- Step 5** Pipette **1.0 ml** of Serum Diluent followed by **10 µl** of Conjugate into each well. Incubate **30 minutes** (±5 min) at room temperature on rocker or rotating shaker.
- Step 6** Repeat **Step 4**.
- Step 7** Pipette **1.0 ml** Substrate into each well and incubate with gentle shaking **5 minutes** (±2 min) at room temperature and reduced light.
- Step 8** Aspirate liquid in wells into waste container. Thoroughly wash strips with deionized water by squirting approximately 2ml of solution directly onto strips. Wash strips with gentle agitation for **5 minutes** aspirate liquid into waste container. **Repeat**.
- Step 9** Using blunt forceps, remove strips from assay tray and place them gently onto absorbent paper. Handle strips only at the ends and let them dry **15-20 minutes**.

Quality Control

Though the control cards are lot specific, negative and positive controls must be included in each test run to ensure proper performance of the assay.

Positive Control Reaction: The positive control strip will exhibit a dense blue-violet band at 30kD, sometimes a doublet (*Figure 1*).

Negative Control Reaction: Typically, the negative control will exhibit no band. Sometimes some sera may give reactions on the strip but not at 30kD.

RESULTS

Reading and Interpretation Guidelines

- Step 1 Hold test strip between the positive and negative control reactions on the provided, laminated Control Card and align test strip using the alignment marker as the reference point.
- Step 2 Compare reaction of test strip with those of the controls on either side. Use of a magnifying glass can facilitate proper alignment of test strips along the molecular weight marker and assist in observation of weak reactions.
- Step 3 Check for a band on the test strip that aligns with the 30kD band on the positive control strip. Such a reaction should be considered positive. Positive reactions can occur in varying intensities from weak to strong. Weak reactions should be compared with baseline reaction intensities at the corresponding position on the negative control strip. Reactions that are distinctly darker or denser than the baseline intensity should be considered positive. Occasionally, minor bands may be visible with certain sera. However, such reactions are not relevant and should be considered as negative for antibodies to P0 (30kD) antigen.

Note: Automated Western blot machines can occasionally show extra bands that may not be related to the 30kD band. Ensure the washing cycle is thorough. Repeat process manually for clarity.

Antibodies to P0 (30kD) antigen occur in patients with autoimmune hearing loss.

Anti-P0 (30kD) Antibodies in Patients with SNHL⁸

Disease Group	No. Tested	No. Positive	% Positive
Sudden	9	1	11
Progressive	20	14	60
Unspecified	5	2	40
Normal	90	1	1

LIMITATIONS OF THE PROCEDURE

The Immublot™ anti-P0 (30kD) antibody Western Blot assay should be used as an aid to the diagnosis of autoimmune hearing loss. Positive results may be found in other autoimmune conditions such as in peripheral neuropathy. Hence results should be evaluated and interpreted by a medical authority in light of the patient's clinical history and other laboratory findings. Some sera may react to antigens of other molecular sizes. These are not of any clinical relevance.

TROUBLESHOOTING GUIDE

- Strong band at 30kD on Negative Control strip. **Likely cause: contaminated Negative Control vial or cross contamination from well containing a positive serum.**
- Positive Control appears like Negative Control strip. **Likely cause: Negative Control vial was confused as Positive Control vial.**
- Strips are completely blank. **Likely cause: addition of Conjugate or Substrate was omitted.**
- High background and poor contrast between bands and background. **Likely cause: wash step (s) may have been omitted or incorrectly performed, or incubations were overextended.**

EN


Figure 1: Diagram of PO Positive Reaction

ImmuBlot™ Anti-P0 Antibody Western Blot

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IT

REFERENCES

1. McCabe BF, Autoimmune sensorineural hearing loss. *Ann Otol.* 1979; 88:5852-589.
2. Arnold W, Pfaltz R, Alternatt HJ. Evidence of serum antibodies against inner ear tissues in the blood of patients with certain sensorineural hearing disorders. *Acta Otolaryngol (Stockh).* 1985; 99:437-444.
3. Harris JP, Sharp PA. Inner ear autoantibodies in patients with rapidly progressive sensorineural hearing loss. *Laryngoscope.* 1990; 100:516-524.
4. Yamanobe S, Harris JP. Inner ear-specific autoantibodies. *Laryngoscope.* 1993; 103:319-325.
5. Shin S, Billings P, Keithley EM, Harris JP. Comparison of anti-heat shock protein 70 (anti-shp-70) and anti-68-kD inner ear protein in the sera of patients with meniere's disease. *Laryngoscope.* 1997; 107:222-227.
6. Harris JP and Aframian D. Role of autoimmunity in contralateral delayed endolymphatic hydrops. *Am J Otol.* 1994; 15:710-716.
7. Yoo TJ, Du X, Kwon SS. Molecular mechanism of autoimmune hearing loss. *Acta Otolaryngol Suppl.* 2002;548:3-9. Review.
8. Tomasi JP, Lona A, Deggouj N, Gersdorff M. Autoimmune sensorineural hearing loss in young patients: an exploratory study. *Laryngoscope.* 2001;11:2050-3.
9. Boulassel MR, Deggouj N, Tomasi JP, Gersdorff M. Inner ear autoantibodies and their targets in patients with autoimmune inner ear diseases. *Acta Otolaryngol.* 2001;121:28-34.
10. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].

For technical assistance please contact:

 **IMMCO Diagnostics, Inc.**
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rev: APR2017

Document No. PI4192