EN

TMB Enzyme Substrate contains an irritant that may be harmful if inhaled, ingested or absorbed through the skin. Do not ingest and avoid contact with skin and eyes.

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

Materials provided

ImmuLisa™ GBM Antibody ELISA **REF** 5154

Kits contains sufficient reagents to perform 96 determinations.

12 x 8	MICROPLATE GBM	Microplate with individual breakaway microwells. Coated with GBM antigen. Ready for use.
1 x 1.75 ml	CONTROL + GBM	Ready to use Positive Control (red cap). Contains human serum positive for GBM antibodies. The expected concentration range in EU/ml is printed on the label.
1 x 1.75 ml	CONTROL -	Ready to use Negative Control (white cap). Contains human serum.
5 x 1.75 ml	CALIBRATOR A GBM CALIBRATOR B GBM CALIBRATOR C GBM CALIBRATOR D GBM CALIBRATOR E GBM	Ready to use set of 5 Calibrators . Calibrator A (green cap) 160 EU/ml, Calibrator B (violet cap) 80 EU/ml, Calibrator C (blue cap) 40 EU/ml, Calibrator D (yellow cap) 20 EU/ml, and Calibrator E (orange cap) 1 EU/ml. Derived from human serum containing GBM antibodies. Concentrations in EU/ml are printed on the labels.
1 x 15 ml	IgG-CONJ HRP	HRP goat anti-human IgG Conjugate . Ready for use. Color coded pink.
1 x 60 ml	DIL	Serum Diluent. Ready for use. Color coded purple.
1 x 15 ml	SUBSTRATE TMB	TMB enzyme substrate. Ready for use. Protect from light.
1 x 15 ml	STOP H2SO4	Stop Solution*. Ready for use.
2 x vials	BUF WASH	Powder Wash Buffer. Reconstitute to one liter each.
1 x		Protocol Sheets

Optional Components

1 x 60 ml **BUF** **WASH** Liquid concentrated Wash Buffer. **Reconstitute to one liter.**

Symbols used on labels:

LOT Lot number

REF Catalog number

IVD In vitro diagnostic use

 Use by

 Storage temperature

 Consult instructions for use

 Number of tests

 Manufacturer

 Date of Manufacture

 *Danger; Causes severe skin burns and eye damage. Causes severe eye damage. Wear protective gloves, protective clothing, eye protection, face protection. Wash exposed skin thoroughly after handling. If in eyes: rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician.

Material required but not provided

- Deionized or distilled water

EN

- Squeeze bottle to hold diluted wash buffer
- Pipettes capable of delivering 5 µl to 1000 µl
- Disposable pipette tips
- Clean test tubes 12 x 75 mm and test tube rack
- Timer
- Absorbent paper towels
- Microplate reader capable of reading absorbance values at 450 nm. If dual wavelength microplate reader is available, the reference filter should be set at 600–650 nm
- Automatic microplate washer capable of dispensing 200 µl

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Store specimens at 2°–8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples. It is recommended that frozen specimens be tested within one year

PROCEDURE

Procedural Notes

- Carefully read the product insert before starting the assay.
- All dilutions of the patient samples should be prepared prior to starting with the assay.
- Let patient specimens and test reagents equilibrate to room temperature before starting with the test procedure. It is suggested that reagents be left on the bench out of the box for 30 minutes prior to use. Return all unused specimens and reagents to refrigerator immediately after use.
- Remove required microwell strips from the pouch and carefully reseal the pouch to prevent condensation in the unused wells. Return pouch immediately to refrigerator.
- **Good washing technique is critical.** If washing is performed manually, adequate washing is accomplished by directing a forceful stream of wash buffer with a wide tip wash bottle across the entire microplate. **An automated microplate washer is recommended.**
- Use a multichannel pipette capable of delivering 8 or 12 wells simultaneously. This speeds the process and provides more uniform incubation times.
- For all steps, careful control of timing is important. The start of all incubation periods begins with the completion of reagent addition.
- Addition of all samples and reagents should be performed at the same rate and in the same sequence.

Test Method

Step 1 Let all reagents and specimens equilibrate to room temperature.

Step 2 Label protocol sheet to indicate sample placement in the wells. It is good laboratory practice to run samples in duplicate.

Step 3 For a **qualitative determination** use Calibrator D (*vial with yellow cap*) only.
or

For a **semi-quantitative determination** use Calibrators A through E as depicted in the sample layout below.

Qualitative			
A	Blank	S5	Etc.
B	-	S6	
C	+	S7	
D	Cal D	S8	
E	S1	S9	
F	S2	S10	
G	S3	S11	
H	S4	S12	
	1	2	3

Semi-Quantitative			
A	Blank	S1	Etc.
B	-	S2	
C	+	S3	
D	Cal A	S4	
E	Cal B	S5	
F	Cal C	S6	
G	Cal D	S7	
H	Cal E	S8	
	1	2	3

Step 4 Prepare a 1:101 dilution of the patient samples by mixing 5 µl of the patient sera with 500 µl of Serum Diluent.

EN

- Step 5** Remove the required microwells from pouch and return unused strips in the sealed pouch to refrigerator. Securely place the microwells into the extra provided holder.
- Step 6** Pipette **100 µl** of Ready to use Calibrators, Positive and Negative controls and diluted patient samples (**1:101**) to the appropriate microwells as per protocol sheet.
Note: Include one well which contains **100 µl** of the Serum Diluent as a reagent blank. Zero the ELISA reader against the reagent blank.
- Step 7** Incubate **30 minutes** (±5 min) at room temperature.
- Step 8** Wash **4x** with wash buffer. For manual washing, fill each microwell with reconstituted wash buffer. Discard the fluid by inverting and tapping out the contents of each well or by aspirating the liquid from each well. To blot at the end of the last wash, invert strips and tap the wells vigorously on absorbent paper towels. For automatic washers, program the washer as per manufacturer's instructions.
- Step 9** Pipette **100 µl** of Conjugate into microwells.
- Step 10** Incubate **30 minutes** (±5 min) at room temperature.
- Step 11** Wash all microwells as in Step 8.
- Step 12** Pipette **100 µl** of Enzyme Substrate into each microwell in the same order and timing as for the Conjugate.
- Step 13** Incubate **30 minutes** (± 5 min) at room temperature.
- Step 14** Pipette **100 µl** of Stop Solution into each microwell using the same order and timing as for the addition of the Enzyme Substrate. Read absorbance values within **30 minutes** of adding Stop Solution.
- Step 15** Read absorbance of each microwell at **450 nm** using a single or at 450/630 nm using a dual wavelength microplate reader against the reagent blank set at zero absorbance.

Quality Control

Calibrators, Positive and Negative Controls and a reagent blank must be run with each assay to verify the integrity and accuracy of the test. The absorbance reading of the reagent blank should be <0.3. The Calibrator A should have an absorbance reading of not less than 1.0, otherwise the test must be repeated. The negative control must be <10 EU/ml. If the test is run in duplicate, the mean of the two readings should be taken for determining EU/ml. While performing Qualitative determinations, the optical density of Calibrator D must be greater than that of the negative control and less than the absorbance of the positive control. For semi-quantitative determinations the positive control must give values in the range stated on the vial.

RESULTS

Calculations

The concentrations of the patient samples can be determined by either of two methods:

1. QUALITATIVE DETERMINATION

$$\frac{\text{Abs. of Test Sample}}{\text{Abs. of Calibrator D}} \times \text{EU/ml of Calibrator D} = \text{EU/ml Test Sample}$$

It is recommended that qualitative results be reported as "positive" or "negative." Sample results greater than or equal to Calibrator D are considered positive.

2. SEMI-QUANTITATIVE DETERMINATION

Plot the absorbance of Calibrator A through E against their respective concentrations on linear-log graph paper. Plot the concentrations in EU/ml on the X-axis against the absorbance on the Y-axis and draw a point-to-point curve fit. Determine the concentrations of the patient samples from the curve according to corresponding absorbance values. Alternately, a four parameter curve may be used to plot the standard curve.

It is recommended that semi-quantitative results be reported as "positive," "negative," or "indeterminate" with EU/ml unit values. Indeterminate/ borderline results should be retested and evaluated along with other laboratory methods.

Interpretation

Interpretation values were determined by testing 64 normal blood donors and disease control specimens. The cutoff was established using mean of the normal subjects plus 3.5 SD for the GBM Ab ELISA and assigned an arbitrary value of 20 EU/ml. IMMCO suggests use of the reference range below. Each laboratory should validate assay values for their own conditions.

GBM Ab value	Interpretation
<20 EU/ml	Negative
20-25 EU/ml	Indeterminate (Borderline)
>25 EU/ml	Positive

Calibrator

The Ready to Use Calibrators are included to provide semi-quantitation and must be used with each run. Patient samples containing high antibody levels may give absorbance values greater than that of Calibrator A. For determining accurate semi-quantitative values, such specimens should be further diluted so they fall within the range of the calibrator curve when retested. For determining EU/ml values, multiply the units obtained by the dilution factor.

LIMITATIONS OF PROCEDURE

The assay should not be performed on grossly hemolyzed, microbially contaminated or lipemic samples. This method should be used for testing human serum samples only. The results obtained serve only as an aid in the diagnosis. Taken alone, these results should not be interpreted as diagnostic. Serum of patients with RPGN may be negative for GBM antibodies. Store specimens at 2–8°C for no longer than one week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

EXPECTED VALUES

Test results in a normal population are usually expected to be negative. Studies of clinical populations were tested to survey positive incidence for GBM antibodies with this assay. Populations included subjects with Goodpasture's syndrome, autoimmune disease controls and presumed disease-free normal human subjects.

Disease	Expected	n	n Pos	% Pos
Goodpasture's Syndrome	95%	65	61	93.8%
Systemic lupus erythematosus	<5	8	0	0.0%
Sjögren's syndrome	<5	16	0	0.0%
Systemic sclerosis	<5	8	0	0.0%
Thyroiditis	<5	16	0	0.0%
Mixed connective tissue disease	<5	8	0	0.0%
Rheumatoid arthritis	<5	8	0	0.0%
Celiac disease	<5	8	0	0.0%
Vasculitis	<5	8	0	0.0%
Normal Human Sera	<5	163	4	2.5%
Total NHS and Disease Controls		243	4	1.6%

PERFORMANCE CHARACTERISTICS

The utility of the Immulisa™ GBM Antibody ELISA was evaluated by testing Goodpasture's specimens alongside disease controls and "normal" human sera. These specimens were also tested on commercially available test kits. Only specimens in the linear range of the assay were included in the method comparison. These results are summarized below.

A. Immulisa™ GBM Antibody ELISA vs. other GBM Antibody kit:

		Other GBM Antibody ELISA		
		Positive	Negative	Total
IMMCO	Positive	42	5	47
GBM Ab	Negative	0	79	79
ELISA	Total	42	84	126

Positive Percent Agreement: 100% (95% CI 89.6%–100%)
 Negative Percent Agreement: 94% (95% CI 86.0%–97.8%)
 Relative Agreement: 96% (95% CI 90.5%–98.5%)

B. Cross Reactivity: A total of 68 from individuals with associated connective tissue disorders or other potentially cross-reactive autoimmune disorders were selected to test for GBM antibodies using the Immulisa™ assay.

Condition	n	Positive
Celiac disease	8	0
Hashimoto's disease	16	0
Rheumatoid arthritis	8	0
Sjögren's syndrome	16	0
Systemic lupus erythematosus	8	0
Vasculitis	8	0
Total	64	0 (0%)

EN Precision

Precision was tested with positive specimens selected throughout the range of the assay. Assay runs of 6 replicates of each specimen were conducted on 13 days. Repeatability was determined with 12 replicates of each specimen.

S #	Mean (EU/ml)	Total Imprecision		Between days		Within run (Repeatability)	
		SD (EU/ml)	CV%	SD (EU/ml)	CV%	SD (EU/ml)	CV%
1	12.0	1.731	14.5%	1.769	14.9%	1.521	12.9%
2	14.9	1.438	9.7%	1.461	10.0%	1.284	8.7%
3	27.0	2.536	9.4%	2.685	9.6%	1.501	5.4%
4	31.4	3.117	9.9%	3.113	10.5%	2.689	9.1%
5	106.5	4.743	4.5%	4.271	4.1%	2.914	2.8%
6	148.7	7.456	5.0%	5.902	4.1%	7.544	5.2%

Reproducibility

Qualitative reproducibility was tested with 90 runs of samples in the negative range, ~20% below cutoff, ~20% above cutoff and in the moderate positive range of the assay using the qualitative analysis method. 6 replicates were tested in 13 runs and an additional 12 replicates were tested in a following run over multiple days. Assay results for these tests produced 100% qualitative (positive/negative) agreement.

Limit of Detection

The limit of detection (LoD) was determined based on 60 replicates of the blank and 10 replicates each of 6 low-level (NHS) samples. LoD was determined to be 2.4 EU/ml.

Linearity and Recovery

Linearity and recovery were tested by diluting positive specimens through the assay range in equidistant dilutions and comparing actual vs. expected results. The linear range of the assay was determined to be 2.4 (LoD) – 160 EU/ml. Results are summarized below:

Test Range (EU/ml)	Slope (95% CI)	Y-intercept (95% CI)	R ²	% recovery (obtained/expected)
1.9 to 92.3	1.01 (.923 to 1.105)	0.7 (-4.2 to 5.6)	0.992	91% to 116%
19.7 to 108.5	1.01 (.891 to 1.131)	2.9 (-5.1 to 10.9)	0.986	101% to 117
31.1 to 180.0	0.98 (.889 to 1.070)	1.2 (-9.5 to 11.8)	0.992	92% to 104%

Interference

Interference was studied by mixing sera with known GBM antibody levels with potentially interfering serum samples and studying deviation from expected results. No significant interference was demonstrated for the following substances at the levels indicated: hemoglobin (2 g/L), bilirubin (342 µmol/L), rheumatoid factor (100 EU/ml), cholesterol (13 mmol/L), and triglycerides (37 mmol/L).

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1. Jayne DRW, Marshall PD, Jones SJ and Lockwood CM. Autoantibodies to GBM and neutrophil cytoplasm in rapidly progressive glomerulonephritis. *Kidney Int*; 1990, 37:965-970.
2. Saxena R, Bygren P, Arvastson B, Wieslander J. Circulating autoantibodies as serological markers in the differential diagnosis of pulmonary renal syndrome. *J Intern Med*. 1995. 238:143-152.
3. Wieslander J, Bygren P, Heinegård D. Anti-basement membrane antibody: Immunoenzymatic assay and specificity of antibodies. *Scand J Clin Lab Invest*. 1981. 41:763-772.
4. Butkowski R, Langeveld J, Wieslander J, Hamilton J, Hudson BG. Localization of the Goodpasture epitope to a novel chain of basement membrane collagen. *J Biol Chem*, 1987; 262:7874-7877.
5. Hellmark T, Johansson C, Wieslander J. Characterization of anti-GBM antibodies involved in Goodpasture's Syndrome. *Kidney Int*. 1994. 46:823-829.
6. Segelmark M, Butkowski R, Wieslander J. Antigen restriction and IgG subclasses among anti-GBM autoantibodies. *Nephrol Dial Transplant*. 1990. 5:991-996.
7. Wieslander J, Bygren P, Heinegård D. Isolation of the specific glomerular basement membrane antigen involved in Goodpasture Syndrome. *Proc Natl Acad Sci USA*. 1984. 81:1544-1548.
8. Wieslander J, Barr JF, Butkowski RJ, Edwards SJ, Bygren P, Heinegård D, and Hudson BG. Goodpasture antigen of the glomerular basement membrane: localization to noncollagenous regions of type IV collagen. *Proc Natl Acad Sci USA*. 1994. 81:3838-3842.
9. Hudson BG, Wieslander J, Wisdom B, Noelken ME. Biology of disease. Goodpasture Syndrome: Molecular architecture and function of the basement membrane antigen. *Lab Invest*. 1989. 61:256-269.
10. Hudson BG, Reeders ST, Tryggvason K. Type IV collagen: Structure, gene organization, and role in human disease. *J Biol Chem*. 1993. 268:26033-26036.
11. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control National Institutes of Health. 2007; (HHS Pub No [CDC] 93-8395)

For technical assistance please contact:

 **IMMCO Diagnostics, Inc.**
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rev. APR2017

Document No. PI5154