ImmcoStripe™ 68 kD (hsp-70) Antibody Line Immunoassay (LIA)

IVD

PRODUCT INSERT

REF 6001 68 kD (hsp-70) Antibody LIA 20 Determinations

INTENDED USE

A line immunoassay for use in the detection and identification of antibodies to the 68 kD (hsp-70) antigen associated with sensorineural hearing loss (SNHL).

SUMMARY AND EXPLANATION

Hearing loss may be caused by a number of conditions. Certain types of hearing loss can be reversed if diagnosed early and appropriate treatment is instituted. Sensorineural hearing loss (SNHL), commonly referred to as nerve deafness, may be caused by genetic or acquired factors such as infections or can be immunologically mediated. Correct diagnosis of SNHL can be made with the aid of a combination of comprehensive patient history and laboratory studies. In the majority of cases, no cause of SNHL is apparent: such cases are referred to as idiopathic SNHL. A subgroup of idiopathic SNHL cases is treatable with immunosuppressive therapy with gratifying results.^{1,2} The laboratory work up to identify these cases should include serum antibody tests to 68 kD (hsp-70) inner ear antigen.³⁻⁸ Twenty-two percent of patients with bilateral rapidly progressive SNHL and 30% of the patients with Meniere's disease had antibodies that reacted with 68 kD antigen present in the bovine inner ear extracts.⁷ Anti-68 kD (hsp-70) antibodies also occur in approximately 60% of patients with bilateral and 35% of patients with unilateral Meniere's Syndrome and 37% of patients with contralateral delayed endolymphatic hydrops. In situations where corticosteroids are contraindicated, methotrexate or cytoxan treatment can be prescribed.⁹

PRINCIPLES OF PROCEDURE

This line immunoassay uses purified recombinant inducible hsp-70 antigen (rhsp-70) from bovine kidney. The hsp-70 antigen is immobilized onto nitrocellulose membrane strips along with conjugate, serum and cut-off control lines. The serum and conjugate control lines serve as internal procedure controls to verify the addition of serum and conjugate, respectively. The cut-off line provides a colorimetric standard for the evaluation of the reacted hsp-70 test line.

To perform the test, strips are incubated with diluted patient serum. In positive sera, antibodies specifically bind to the rhsp-70 protein on the strip. The strips are washed according to the protocol and the reconstituted conjugate is added to the test strips. After incubation and wash steps, the ready to use substrate is added to the strips. During a 10 minute incubation, conjugate and substrate binding produces visible violet lines for serum, conjugate and cut-off control lines. If the sample is positive for anti-hsp-70 antibodies, the test line will show a reaction more intense than the cut-off line. Reactions are read visually and reported as positive, negative or equivocal (comparable to cut-off line).

Materials Provided

68 kD (hsp-70) Antibody LIA

REF 6001

Kits contains sufficient reagents to perform 20 determinations.

20 x STRIP|HSP70|LIA

Line Immunoassay Test Strips, containing antinuclear antigen coated test lines and control lines. Ready for use.

1 x 120 µl CONTROL|+|HSP70

Positive Control (red cap). Contains human serum positive for antibodies against HSp70 antigen.

1 x 30 ml CONJ|LIA

IgG Conjugate.

1 x 30 ml SUBSTRATE

Enzyme Substrate (amber bottle). Ready for use. **Protect from light**.

1 x 50 ml SAMPLE|DIL

Sample Diluent

1 x 50 ml BUF|WASH|LIA

Wash Buffer Concentrate. Reconstitute to one liter with deionized or distilled water or as needed proportionally.

2 x LIA 10 well Assay Trays

Symbols used on labels:

	Lot number
	Catalog number
	In vitro diagnostic use
	Use by
	Storage temperature
	Consult instructions for use
	Number of tests
	Manufacturer
	Date of Manufacture

Materials Required But Not Provided

- Clean 1.000 ml graduated cylinder
- Non-serrated forceps (Filter forceps)
- Rocker or rotating platform shaker
- Absorbent paper or paper towels
- Deionized or distilled water
- Squeeze bottles to hold diluted wash buffer or distilled water
- Pipettes capable of delivering 10 to 1.000 µl
- Disposable pipette tips
- Timer

REAGENTS**Materials provided**

Store all reagents at 2–8°C; **do not freeze.**

All reagents must be brought to room temperature (18–25°C) and mixed thoroughly prior to use. Do not use if reagent is not clear or if insoluble precipitate is present. After reconstitution, the reagents are stable until the indicated expiry date when stored at 2–8°C and protected from contamination, unless state otherwise below:

- Antigen coated test strips **[STRIP|HSP70|LIA]** are ready for use. Please allow the test strip bag to reach the room temperature before opening to avoid condensation and associated deterioration. Please re-pack unused test strips and store at 2–8°C in dark and dry conditions.
- **Sample Diluent** is ready to use **[SAMPLE|DIL]** is stable for at least 8 weeks when stored properly and protected from microbial or chemical contamination.
- Reconstitute 1 part **[BUF|WASH|LIA]** into 19 parts of distilled or deionized water to result in 1 liter of **Wash Buffer**. Wash Buffer is stable for at least 8 weeks after reconstitution when stored properly and protected from microbial contamination.
- **[CONJ|LIA]** and **[SUBSTRATE]** are stable for at least 8 weeks after opening when stored properly and protected from microbial contamination. **[SUBSTRATE]** is light sensitive and must be stored in the provided amber colored bottle.

Antigen strips can only be used once. Do not interchange components of different lots. Do not use reagents beyond expiration date indicated on labels.

Precautions

All human derived components used have been tested for HBsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. However, human blood derivatives and patient specimens should be considered potentially infectious. Follow good laboratory practices in storing, dispensing and disposing of the above materials.¹⁰

WARNING – Proclin 300 is a preservative. Upon disposal of liquids containing Proclin 300, flush with large volumes of water to dilute the components below active levels.

EN

Instructions should be followed exactly as they appear in this kit insert to ensure valid results. Do not interchange kit components with those from other sources. Follow good laboratory practices to minimize microbial and cross contamination of reagents when handling. Do not use kit components beyond expiration date on the labels.

SPECIMEN COLLECTION AND HANDLING

Only serum specimens should be used in this procedure. Specimens with gross hemolysis, elevated lipids or microbial contamination may interfere with the performance of the test and therefore must not be used. Store specimens at 2–8°C for no longer than one week. For longer storage, serum specimens should be frozen. It is recommended that frozen specimens be tested within one year. Avoid repeated freezing and thawing of samples.

PROCEDURE

Procedural Notes

- Read Product Insert carefully before starting with the assay.
- Let serum specimens and test reagents equilibrate to room temperature for approximately 30 minutes prior to starting the test procedure. Return all unused specimens and reagents to the refrigerator promptly after use.
- Proper washing technique is critical to the satisfactory performance of the assay.
- Manipulate test strips with clean forceps or gloves only. Avoid touching the white nitrocellulose areas.
- The test lines are placed above (refer to Figure 1 Schematic) the cut-off, serum and conjugate control lines as described in the schematic. Serum and conjugate control lines appear on the same piece of nitrocellulose at the bottom.
- Assign specimen identification numbers to the respective strips on the Report Form. Each strip has strip number and lot number printed on the bottom for traceability.
- Complete all other relevant information on the Report Form prior to starting the assay.

Test Method

- Step 1** Using gloves or blunt forceps, peel off the required number of strips. Care should be taken not to touch the nitrocellulose coated areas with bare hands or pointed forceps.
- Step 2** Place required number of strips **STRIP|HSP70|LIA** labeled side up into individual wells of the assay tray.
- Step 3** Pipet 1.5 ml of Sample Diluent **SAMPLE|DIL** into each well making sure that the strips are completely submerged under the liquid.
- Step 4** Incubate the strips in Sample Diluent **SAMPLE|DIL** for at least 10 minutes. The blue tracking line starts to disappear as the membrane is soaked.
- Step 5** Pipet 15 µl of serum or control sample into appropriate wells to obtain a 1:101 dilution. Incubate 60 minutes (±5 min.) at room temperature on a rocker or rotating shaker.
- Step 6** WASH: Aspirate sample solution into waste container. Thoroughly wash strips with Wash Buffer **BUF|WASH|LIA** by squirting approximately 2 ml of solution directly onto strips. Wash strips with gentle agitation for 5 minutes and aspirate solution into waste container. Repeat the wash two more times. Caution: Complete washing of the strips between incubations is crucial to obtain valid results. Improper washing will result in high background staining. Do not allow the strips to become dry at any step during the assay.
- Step 7** Pipet 1.0 ml of conjugate **CONJ|LIA** into each well. Incubate 30 minutes (±5 min) at room temperature on rocker or rotating shaker.
- Step 8** Repeat Step 6.
- Step 9** Pipet 1.0 ml **SUBSTRATE** into each well and incubate with gentle shaking 10 minutes at room temperature and reduced light. The serum and conjugate control lines develop intense color after incubation in substrate. Cut-off control line develops into a negative to faintly colored line after the incubation.
- Step 10** To stop the reaction, rinse strips 2x with distilled/deionized water by squirting approximately 2 ml of water directly onto strips followed by aspiration. Do not soak/wash for more than 10 minutes as this may result in decreased sensitivity of the developed colored lines.
- Step 11** Using blunt forceps remove strips from assay tray and place them gently onto absorbent paper and allow them to dry. Let the strips dry before analysis or affixing them on the report/scoring sheet.

EN Quality Control

Procedural Controls: Each strip has three procedural controls for the addition of serum and conjugate and a cut-off line for determining the weak or negative reactions.

Positive and Negative controls are available as optional components and may be run for additional quality control. Individual labs are expected to optimize the substrate development time by +/- 2 minutes based on the blot processor or methods that they use. It is suggested that the cut-off line should be faintly visible line post incubations with substrate. An invisible or too a strong cut-off line suggests optimization of the substrate development time for the given run conditions.

Interpretation

The test strips contain control lines at the bottom and test line above the controls. The bottom end of the test strip (near serial number) has three control lines: the cut-off line, the serum control line and the conjugate control line from top to bottom. The cut-off allows the technician to determine the test result as positive, negative or indeterminate (+/-). The two procedure control lines ensure the addition of specimen, conjugate and substrate.

Compare the reaction of test line with those of the controls. Use of a magnifying glass can assist in observation of weak reactions.

- As labeled in Figure 1, the serum and conjugate control lines should be clearly positive indicating a successful experiment. The cut-off is a faint line with variation in intensity based on the experimental conditions. The schematic in Figure 1 shows an example test line and 3 control lines. The test line development depends on the sample. Positive reactions can occur in varying intensities from weak to strong. **Weak reactions should be compared with intensity of the provided cut-off line within the strip.** Reactions that are distinctly darker or denser than the intensity of the cut-off line should be considered positive.
- Strips may show a homogeneous or discolored background due to various interfering factors in lipemic or hemolytic sera. This effect can also be seen if the test strips are not sufficiently blocked or accidentally allowed to dry up during the assay.
- In case of weak positive and negative reactions, the reacted line intensity should be compared to cut-off line to determine the result as negative (weaker intensity than the cut-off line) or equivocal (+/-). Indistinguishable from cut-off line).
- Dried strips can be assembled in the provided report/scoring sheet. The plastic protective flap is permanently affixed to the report sheet on the left edge. Carefully peel the plastic flap in the right to left direction like a page of the book. Place the reacted strips on the adhesive tape in the respective slot and cover the plastic flap back in place. The protective plastic flap is designed to be reusable for multiple sessions of experiments and the strips can be assembled in respective slots. The technician can use the form to record the lot numbers of used reagents, specimen number and results/comments.


Figure 1: Schematic of reacted LIA strips with one test line.

Figure 2: Schematic of Report Sheet

LIMITATIONS OF PROCEDURE

Only serum specimens should be used in this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of the test and should not be used. Do not store specimen at 2–8°C more than a week. For longer storage, serum specimens should be frozen. Avoid repeated freezing and thawing of samples.

The 68 kD (hsp-70) Antibody Line Immunoassay should be used as an aid to diagnosis. Positive results may be found in other autoimmune conditions or certain infectious diseases. Hence results should be evaluated and interpreted by a medical authority in light of the patient's clinical history and other laboratory findings. Some sera may react to the MW marker occasionally, the significance of which is not known.

EXPECTED VALUES

Antibodies to 68 kD (hsp-70) occur in patients with active idiopathic SNHL and antibody titers are shown to correlate with disease activity.

EN**Anti-68 kD (hsp-70) Antibodies in Patients with Idiopathic Bilateral SNHL¹¹**

Disease Group	No. Tested	No. Positive	% Positive
IPBSNHL	72	42	58
Otosclerosis	11	0	0
Cogan's Syndrome	8	0	0
Normals	53	1	2

Correlation of Anti-68 kD (hsp-70) Antibody Reactivity to Disease Activity¹¹

Disease Activity	Anti-68 kD Antibody
Active	89%
Inactive	0%

One third-party study of 34 patients with rapidly progressive hearing loss demonstrated that anti-rhsp-70 OTOBlot had 42% sensitivity, 90% specificity, and 91% positive predictive value for predicting steroid responsiveness.¹²

PERFORMANCE CHARACTERISTICS

The utility of the ImmuStripe™ 68 kD (hsp-70) Antibody Line Immunoassay was evaluated by testing well-characterized specimens from SNHL patients, idiopathic hearing loss patients, disease controls and “normal” human sera. These results are summarized below.

		Clinical Status		
		Positive	Negative	Total
68 kD	Positive	36	2	38
(hsp-70)	Negative*	14	79	93
LIA	Total	50	81	131

Sensitivity: 72.0%

Specificity: 97.5%

Clinical Agreement: 87.8%

* Indeterminate samples were considered negative.

SNHL/Idiopathic Hearing Loss Subjects: 50

Disease Controls: 33

Healthy Normal Subjects: 48

Reproducibility

Assays of samples in the negative range, equivocal and positive range were performed to determine qualitative reproducibility from run to run and operator to operator. Results produced 100% qualitative agreement.

Cross Reactivity

Autoimmune disease controls were reacted to test potential cross-reactivity of the assays, including the following patient populations: celiac disease (5), extractable nuclear antigen antibody positive connective disease (12), Hashimoto's thyroiditis (6), rheumatoid arthritis (5), and systemic lupus erythematosus (5). One rheumatoid arthritis specimen was positive yielding a 3.3% positive rate in this population.

Interference

Interference was studied by mixing sera with known 68 kD (hsp-70) antibody results with potentially interfering serum samples and studying deviation from expected results. No significant interference impacting qualitative results was demonstrated for the following substances at the levels indicated: hemoglobin (2 g/L), bilirubin (342 µmol/L), triglycerides (25 mg/ml) and Rheumatoid Factor (100 EU/ml).

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REFERENCES

1. McCabe BF. Autoimmune sensorineural hearing loss. *Ann Otol.* 1979; 88:5852-589.
2. Arnold W, Pfaltz R, Alternatt HJ. Evidence of serum antibodies against inner ear tissues in the blood of patients with certain sensorineural hearing disorders. *Acta Otolaryngol (Stockh).* 1985; 99:437-444.
3. Harris JP, Sharp PA. Inner ear autoantibodies in patients with rapidly progressive sensorineural hearing loss. *Laryngoscope.* 1990; 100:516-524.
4. Yamanobe S, Harris JP. Inner ear-specific autoantibodies. *Laryngoscope.* 1993; 103:319-325.
5. Shin S, Billings P, Keithley EM, Harris JP. Comparison of anti-heat shock protein 70 (anti-shp-70) and anti-68-kD inner ear protein in the sera of patients with meniere's disease. *Laryngoscope.* 1997; 107:222-227.
6. Harris JP and Aframian D. Role of autoimmunity in contralateral delayed endolymphatic hydrops. *Am J Otol.* 1994; 15:710-716.
7. Gottschlich S, Billings P, Keithley EM et al. Assessment of serum antibodies in patients with rapidly progressive sensorineural hearing loss and meniere's disease. *Laryngoscope.* 1995; 105:1347-1352.
8. Billings PB, Keithley EM, Harris JP. Evidence linking the 69 kilodalton antigen identified in progressive sensorineural hearing loss patient sera with heat shock protein 70. *Ann Otol Rhinol Laryngol.* 1995; 104:181-188.
9. Sismanis A, Thompson T, Willis HE. Methotrexate therapy for autoimmune hearing loss: a preliminary report. *Laryngoscope.* 1994; 104:932-934.
10. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].
11. Moscicki RA, San Martin JE, Quintero CH et al. Serum antibody to inner ear proteins in patients with progressive hearing loss – correlation with disease activity and response to corticosteroid treatment. *JAMA.* 1994; 272:611-616.
12. Hirose K, Werner MH, and Duckert LG. Utility of laboratory testing in autoimmune inner ear disease. *Laryngoscope.* 1999; 109:1749-1753.

For technical assistance please contact:

immco DIAGNOSTICS, Inc.
60 Pineview Drive
Buffalo, NY 14228-2120 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: (800) 537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands